# Carisolv™ power drive

DE	Bedienungsanleitung	3
GB	Instructions for use	11
IT	Istruzioni per l'uso	19
NL	Gebruiksaanwijzing	27
PL	Instrukcja Stosowania	35
SE	Bruksanvisning	43

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## Contents:

1	INTRODUCTION .....	12
	1.1 General information .....	12
	1.2 Components included .....	12
	1.3 Control unit .....	12
	1.4 Tips .....	13
	1.5 Motor unit with contra-angle .....	13
2	CLINICAL OPERATION .....	13
	2.1 Preparations .....	13
	2.2 Use .....	13
	2.3 Disinfection, cleaning and sterilisation .....	14
3	TECHNICAL OPERATION, CONNECTION .....	15
	3.1 Motor unit with contra-angle .....	15
	3.2 Control unit .....	16
	3.3 Safety/waste and storage .....	17
	3.4 Maintenance .....	17
4	TECHNICAL SPECIFICATIONS .....	18
5	EXPLANATION OF SYMBOLS .....	50
6	MANUFACTURER .....	51

Please read these instructions very carefully.  
You and all employees should be familiar with the product before it is used.

# 1 INTRODUCTION

## 1.1 General information

The Carisolv™ power drive is designed to be used together with Carisolv™ gel and the tips are specially designed for Carisolv™ power drive. Using this tool without Carisolv™ gel is not recommended. The Carisolv™ power drive can be used in combination with Carisolv™ hand instruments.

Carisolv™ is a modern method of treatment for patient-friendly and tissue-preserving removal of caries. The method's main areas of application are for treating root caries, deep cavities, for children or nervous patients, and in all cases where anaesthetics cannot be given.

The method has two primary advantages:

- Considerably less pain and discomfort for the patient
- It helps the dentist to preserve the maximum amount of healthy dentine

The main advantages of the Carisolv™ power drive are:

- Fixed, low rotation speed, which minimises noise level and discomfort
- Limited torque, which helps the dentist to protect healthy dentine
- Specially designed tips help the dentist to work both efficiently and in a way that preserves tissue
- Clear light indicators show the selected mode, minimising the risk of operating errors
- Portable – the device can be moved from room to room or to other clinics

## 1.2 Components included (see the picture, page 2)

1. Motor unit with contra-angle	7. Mounting bracket
2. Control unit with integrated rechargeable batteries	8. Spray adapter
3. Battery charger	9. Protective sleeve for the motor
4. 2 tips Ø 1.0 mm, length 21 mm (not in the picture)	10. Operating instructions (not in the picture)
5. 2 tips Ø 1.5 mm, length 21 mm (not in the picture)	11. Foot control (not included in the kit, available separately)
6. 2 tips Ø 2.0 mm, length 21 mm (not in the picture)	

## 1.3 Control unit

The control unit contains rechargeable batteries, which for safety reasons cannot be charged while they are in use. The motor, foot control and battery charger can be connected to the control unit. The motor, which drives the contra-angle, is controlled by a microprocessor inside the control unit.

*Control unit*

## 1.4 Tips

The tips, or working parts, are specially designed. They are designed to work in two directions; clockwise and anticlockwise. When rotating clockwise (red mode), the cutting edges are more active than when the tip rotates anticlockwise (green mode). These characteristics of the tips, coupled with the fact that the processor in the control unit measures the torque used, help the dentist to both efficiently dissolve and remove the caries with gentle scraping. There is a green mark on the handle of the tips.

*The Carisolv™ power drive tips*

## 1.5 Motor unit with contra-angle

The motor and the contra-angle are built together as one unit, where the motor can be pulled rearward to enable cleaning and autoclaving of the contra-angle.

# 2 CLINICAL OPERATION

## 2.1 Preparations

The Carisolv™ power drive is assembled according to the instructions in “Technical operation”, chapter 3.

If necessary, the cavity should be opened with a conventional rotating instrument or a manual instrument. Any remaining filling that is not to be saved, should be removed.

## 2.2 Use

- 2.2.1 Carisolv™ gel should be mixed and the desired amount of gel should be placed in a suitable cup. For complete information about Carisolv™ gel, please refer to the instructions and product information enclosed in the gel package.
- 2.2.2 Select a tip according to the size of the cavity. Using the tip (or a Carisolv™ hand instrument) to apply the gel to the decayed tooth substance. Ensure that the decayed area is soaked with gel.
- 2.2.3 Allow the gel to work for approximately 30 seconds.
- 2.2.4 Begin removing the caries by working clockwise with the more aggressive-cutting work mode (red indicator) to remove the softened surface dentine. Keep the tip tilted/angled against the caries to achieve the fastest possible removal of the softened caries. **Do not press too hard.** If the pressure on the Carisolv™ power drive becomes too high, the drive will stop to limit the risk of removing healthy dentine.

**NOTE!** When the Carisolv™ power drive rotates clockwise – red indicator – the removal process is more intensive and is best suited when starting the caries removal and for more mineralised caries.

When working near healthy dentine or close to the pulp, switch to anticlockwise rotation – green indicator. The tips will then only scrape, which makes it easier for the dentist to preserve healthy dentine.

When using the anticlockwise mode (green), be sure to be very gentle and do not use too much pressure. In the event of excessive pressure or if the treatment is performed too aggressively, the Carisolv™ power drive will stop. This is indicated by red/green flashing. When the pressure is lightened up, the Carisolv™ power drive will automatically start to work again.

- 2.2.5 Remove the decayed dentine and apply more gel gradually and continue removing the softened caries. **Be sure that there is an ample amount of gel in the cavity at all times.** Avoid drying or rinsing the cavity. Switch to anticlockwise rotation (green indicator) when working near healthy dentine and/or near the pulp.  
The Carisolv™ power drive does not need to be switched off when gel is removed or applied to the cavity. Please note that the Carisolv™ power drive sometimes needs to be complemented by the special hand instruments that have been developed for use with Carisolv™ gel.
- 2.2.6 When the cavity appears to be free of caries, examine with a probe without removing the gel. Check particularly for caries at the dentinoenamel junction.
- 2.2.7 If the cavity appears to be free of caries, remove the gel and clean the cavity with a moist cotton pellet. Examine and probe with a sharp probe.
- 2.2.8 If the cavity is not free of caries, apply more gel and continue the treatment.
- 2.2.9 The periphery of the cavity should be adjusted with a suitable instrument. Restore the cavity with the dental filling material of your choice, according to the manufacturer's instructions.

*NOTE! After finished treatment, handle the gel according to the instructions enclosed in the gel package.*

## 2.3 Disinfection, cleaning and sterilisation

- 2.3.1 The contra-angle and the tips** must be cleaned, disinfected and then steam-sterilised at 135°C before being used. **Control unit, motor with cable, protective sleeve and spray adapter cannot be autoclaved.** These components can be disinfected and cleaned by wiping with an alcohol-cleaner impregnated cloth.

The parts of the handpiece and those of the whole unit should not come into contact with strong cleaning agents.

**When wiping the motor or the contra-angle, great care must be taken to ensure that no dirt comes into the gear. The motor axle should always be protected with the protective sleeve when the contra-angle is removed. Never dip the contra-angle into cleaning agents.**

The methods of disinfection and sterilisation must be suitable, verified and validated. The effectiveness against viruses must be proved. The process must correspond to the rules in force and the procedure must be admitted by the appropriated authorities.

### 2.3.2 Contra-angle

Before the contra-angle is **steam-sterilised**, it has to be disinfected and cleaned. It might be helpful to use a soft brush. To protect the contra-angle, spray it with oil spray before the sterilisation.

*Please note that the contra-angle should be lubricated directly after being cleaned. Use KaVo Oil Spray (no. 4119630), W&H Service Oil MD-400 (no.10940001) or an equivalent oil spray intended for dental handpieces. Remove the normal spray nozzle on the spray can and insert the adapter for oil-spray (see the illustration on page 2, point 8) according to the photo on page 15. Remove the contra-angle from the motor according to the instructions in section 3.1 and push the contra-angle over the spray-adapter as shown in the photo. Spray the contra-angle once, for approximately one second, then keep it in the upright position for*

approximately 15 seconds to let the excess oil out into the rim of the adapter. After this it is preferred to wipe off remaining oil at the outside with a cloth. Make sure that no cloth remnants are left inside the housing. **While spraying the contra-angle hot objects and flames should be kept away!** After that, the contra-angle is sealed aseptically into sterile bags and steam-sterilized at 135°C. Immediately after sterilisation, all parts must be removed from the sterilizer. The sterilizer should be checked at regular intervals to ensure that it is in good working order. The components that have been steam-sterilized must be assembled under sterile conditions. Wipe off any overflowing oil.

How to use the spray adapter

### 2.3.3 Tips

The tips are made of a stainless material, commonly used in dental instruments and can be cleaned, disinfected and sterilised according to the procedures used at the clinic for similar equipment.

## 3 TECHNICAL OPERATION, CONNECTION

### 3.1 Motor unit with contra-angle

#### 3.1.1 Disassembly

Unscrew the light-cap of the motor until the threads are completely unscrewed. Pull the motor rearward. **NOTE!** The free motor axle can easily be damaged. It is therefore important that the protective sleeve (see the illustration on page 2, point 9) is fitted over the motor axle when the motor and the contra-angle are not connected.

#### 3.1.2 Cleaning

See above (section 2.3). Make sure that the inner parts of the contra-angle are kept free from Carisolv™ gel and dirt.

#### 3.1.3 Assembly

After the disinfection, cleaning and sterilisation of the contra-angle and the disinfection of the motor, both parts have to be assembled under clean conditions according to normal clinical practice. Pull back the moveable light cap of the motor as far as possible in the direction of the cable. Push the motor in the contra-angle by holding it at the black cable-socket. Then turn the motor in the contra-angle for a few degrees, which helps the crown of the motor to adapt inside the contra-angle. Now screw the cap into the thread of the contra-angle. Note that the cap can be tightened fully. If this cannot be done, it is due to the fact that the motor is not properly in place. Loosen the socket several turns and tighten it again.

#### 3.1.4 Fitting the tips

The tips are fitted after the locking clip has been moved to the side. After fitting the tip, move the locking clip back to locked position. **Always check that the lock is engaged by pulling the tip.**

## 3.2 Control unit

### 3.2.1 Connecting

The connectors for the motor, foot control and battery charger are colour-coded. The battery charger is plugged into the same socket as the motor (green connector). That makes it impossible to use the device when it is connected to the mains.

### 3.2.2 Controls

Green mode: Anticlockwise rotation

Red mode: Clockwise rotation

Stop mode: The machine stops

The panel on the control unit always indicates the selected mode with a red or green light. The motor unit uses the same colours to indicate setting. Furthermore, the motor unit indicates overload with flashing light.

The rotation, clockwise or anticlockwise, are activated by pressing the respective buttons.

### 3.2.3 Foot control

The Carisolv™ power drive can be used with or without the foot control, which is not included in the kit but is available separately. Indicators as above. If the foot control is connected, a short press will change the working mode and the machine starts when the foot control is pressed down.

### 3.2.4 Torque limiter

The Carisolv™ power drive continuously measures torque and automatically stops when the limit is exceeded. The overload is indicated at the motor by flashing red / green light. The Carisolv™ power drive starts automatically again when the load is reduced.

### 3.2.5 Charging

The yellow light on the panel indicates that the batteries are in need of charging. When the batteries are being charged the yellow light flashes and stops when the batteries are fully charged.

For charging connect the battery charger to the green connector and plug it into the mains. Charging stops automatically when the integrated batteries are fully charged. Charging time is approximately 2-3 hours. The Carisolv™ power drive should not be charged if the temperature is above 35°C.

*NOTE! To ensure the longest service life for the integrated batteries, charging should only be done when the yellow indicator light comes on. Also remember that the batteries always must be charged before they become completely uncharged.*

### 3.2.6 Switching off

The unit does not need to be switched off since this is done automatically after a few minutes of inactivity. The unit is reactivated by pressing any of the buttons or the foot control.

### 3.2.7 Mounting bracket

The control unit should ideally be hung in the safety bracket during operation. To make this possible, the safety mounting bracket should be attached upright to a vertical surface; screws are supplied for this purpose.

### 3.3 Safety/waste and storage

To protect against contamination, use sterile, disposable gloves, spectacles, mouth and nose protection, protective clothing and a protective covering for surfaces. Also observe the accident prevention regulations. Used instruments and waste must be stored and transferred in such a way that patients, personnel and the environment cannot be contaminated.

The cleaned, sterile contra-angle of the handpiece and the tools must be stored under suitable conditions. The transport packaging is only intended for transport and non-sterile storage. If the Carisolv™ power drive is not used for a period of time, the batteries should be recharged before operation. When storing the Carisolv™ power drive, take care that all parts and the instructions are kept together. **Before a longer period of storage, the control unit must be charged.** The device must be stored in a temperature range of between 10°C and 40°C. A high level of air humidity must be avoided.

The control unit contains batteries. Never throw the unit into a fire or water. The usage of Carisolv™ power drive in rooms with an explosive gas atmosphere, such as operating theatres, is not permitted. The Carisolv™ power drive should only be used at room temperatures of between 18°C and 28°C. A high level of air humidity is to be avoided; the device is intended for short-term operation.

### 3.4 Maintenance

Aside from the steps described in section 2.3 above; disinfection, cleaning and sterilisation, the Carisolv™ power drive is maintenance free. The equipment is tested for safety and reliability during the manufacturing process.

If problems or questions arise, contact the supplier or MediTeam Dental AB. The Carisolv™ power drive contains no parts that can be repaired at the clinic. If the device is opened (not the removing of the contra-angle from the motor) or if it is modified, the CE marking will no longer be valid. The warranty may also be expired.

## 4 TECHNICAL SPECIFICATIONS

Model	Carisolv™ 3.01, with handset type BF'(IEC)
RPM	350 r/min (anticlockwise), 450 r/min (clockwise)
Torque	Up to approx. 12 mNm
Settings	Anticlockwise or clockwise. Selected via the control unit or foot control
Programme settings	Two; torque/speed combinations
Control	RxI compensation and torque limiter
Power supply	Two integrated, rechargeable NiMH batteries à 1.2 V, 1600 mAh
Sterilisation	The contra-angle and the tips can be steam sterilised at 135°C. The motor and other components can be wiped clean with disinfectant.
Output voltage	290 mA/5.5 V
Dimensions, control unit	145 x 63 x 35 mm
Dimensions, motor unit	Diameter 20 mm, length approx. 128 mm
Dimensions, tips	Diameter 1.0 mm, 1.5 mm and 2.0 mm. Length 21 mm
Weight, control unit	180 g
Weight, motor unit	64 g incl. power cord
Input voltage, battery charger	90-240V, 50-60 Hz
Output voltage, battery charger	5V

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Carisolv™ is a trademark of MediTeam Dental AB.
- 6 HERSTELLER: Hergestellt durch Schlumbohm OHG im Auftrag von MediTeam Dental AB.
- 6 MANUFACTURER: Manufactured by Schlumbohm OHG for MediTeam Dental AB.
- 6 FABRICANTE: Fabbricato dalla Schlumbohm OHG per la MediTeam Dental AB.
- 6 PRODUCNET: Geproduceerd door Schlumbohm OHG voor MediTeam Dental AB.
- 6 PRODUCENT: Wyprodukowane przez Schlumbohm OHG dla MediTeam Dental AB.